Gamma Ray Spectroscopy

James Amidei

University of Colorado Denver

October 6, 2023

Overview

- Introduction
 Types of Nuclear Decay
 Gamma-Electron Interactions
- ② Detector Scintillator Photomultiplier
- 3 Data Analysis Identifying Peaks Finding Unknowns and Error Propagation
- 4 Special Calculation Finding the Energy of the Compton Edge

Introduction

- Gamma rays are emitted by excited nuclei
- Very high energy
 - Can be as much as 10⁷ times greater than photons released from electron transitions.
- Effectively "invisible" to detector
 - We have to rely on various interactions within absorbing material to gain any insight.

Types of Nuclear Decay

- Alpha decay
 - ${}^{A}_{Z}X \rightarrow {}^{A-4}_{Z-2}X' + {}^{4}_{2}\alpha$
- Beta decay
 - Positive ${}_Z^A X \rightarrow {}_{Z-1}^A X^{'} + {}_1^0 \beta + \nu_e$
 - Negative ${}_Z^A X o {}_{Z+1}^A X^{'} + {}_{-1}^0 eta + ar{
 u}_e$
- Gamma emission
 - ${}^{A}_{Z}X \xrightarrow{\text{relaxation}} {}^{A}_{Z}X' + \gamma$

Gamma-Electron Interactions

- Photo-electric effect
 - Lower energies
- Compton scattering
 - Most common interaction within the energy spectrum of this experiment.
- Pair production
 - At higher energies;
 i.e., ~ 1022 keV, or
 twice the rest energy
 of an electron.




Figure: The relative importance of the three major gamma ray interactions. [Knoll, 2010]

Detector

- Made of sodium iodide (Nal)
- Gamma rays enter and interact with atoms in the crystal lattice
- Gamma energy is transferred to the electrons orbiting each atom, causing them to jump to excited states.
- Electrons de-excite and release visible light photons.
- Visible light photons are captured by photocathode at the entrance to a photomultiplier tube (PMT).

- Initial visible light photon knocks a primary electron from the photocathode into the PMT.
- Primary electron is pulled into an electric field emitted from an anode towards a dynode.
- Electron knocks several additional electrons (roughly 6) out of place upon contact with the dynode.
- The loose electrons are pulled into another, stronger electric field towards another dynode where the process repeats.
- This process occurs ten times in the PMT, which results in $\sim 6^{10}$ electrons being released.




Figure: Diagram of sodium iodide scintillation detector with PMT. [Boeglin, 2023]

- Compton backscatter
 - Compton scatters into detector.
- Compton edge
 - Compton scatters inside detector.
- Full energy peak
 - Full gamma energy absorbed




Figure: Typical spectrum for single gamma ray. [Bailey, 2002]

- Analytically identified the most likely locations of the full energy peaks by attempting to identify smaller, ancillary peaks in the data for known substances.
- Used matplotlib.pyplot, numpy, scipy.optimize and scipy.stats python libraries to write code to return specific values associated with each peak.
- For peaks that were more difficult to identify, we took the ratio of a known channel number and corresponding energy to algebraically solve for a likely channel number and used that estimate to focus our code to specific ranges.
- We took the calculated channel numbers and corresponding energies and used linregress() function from scipy to make a line of best fit.




Figure: Spectrum for Mn-54.






Figure: Peak for single 835 keV gamma ray.



Finding Unknowns and Error Propagation

- Used similar methods to find the peaks for the unknown compound.
- Error for each channel was calculated using the curve_fit function from the scipy library.
- Used the equation for line-of-best-fit to solve for the unknown gamma peak energies.
 - C = mE + b
 - $E = \frac{C b}{m}$
- Found error in the calculated energies using the below equation for the propagation of error.
 - $\sigma_E = \sum_{i=1}^{n=3} \left(\sigma_i \frac{\partial E}{\partial x_i} \right)^{1/2}$
 - Where x_1 , x_2 , and x_3 correspond to m, C, and b respectively.

Finding Unknowns and Error Propagation



Finding Unknowns and Error Propagation

```
Peak 1 is at (E1, C1) = (0.680288864768523 ± 0.028591785162186863, 1037.26734817 ± 0.12423879)

Peak 2 is at (E2, C2) = (1.0835812776221945 ± 0.03783534499638608, 1606.17412676 ± 1.43986508)
```

Figure: Energies and Channel Number of the two unknown gamma rays.

https://atom.kaeri.re.kr/old/gamrays.html

Compton Edge Energy

- Being able to find the energy of the Compton edge could help to more accurately find peaks.
- The energy of the recoiled electron due to Compton scattering can be found using the below equation.

•
$$E_{e^{-}} = h\nu - h\nu^{'} = h\nu \left[\frac{\frac{h\nu}{m_0c^2}(1-\cos\theta)}{1+\frac{h\nu}{m_0c^2}(1-\cos\theta)} \right]$$
 [Knoll, 2010]

- There are two extreme cases we can see from this equation:
 - ① $\theta\cong 0$, in which $h
 u\cong h
 u'$ and the Compton electron has very little energy while the scattered gamma ray retains nearly full energy.
 - 2) $\theta=\pi$, in which the gamma ray backscatters towards its direction of origin and the electron recoils in the direction of incidence with the maximum amount of energy that can be transferred in a single Compton interaction.
- In the case where $\theta=\pi$, this gives us the energy of the Compton edge.

References



Glenn Knoll (2010)

Radiation detection and measurement.

John Wiley and Sons, 2010



Werner U. Boeglin (2023)

Modern Lab Experiments, "22. Scintillation Detector."

Werner U. Boeglin, FIU, 2023



David Bailey (2002)

The Germanium Spectrometer.

David Bailey, 2002